**All-India Institute of Medical Sciences** 

Ansari Nagar, New Delhi-29

(RESEARCH SECTION)

Ref. No. 35/Prop/Anaesth/MK/2019-20/RS

Dated: 04.01.2020

Subject: Purchase of Airtraq SP Video Laryngoscope for the Deptt. of Anaesthesiology,

AIIMS, New Delhi-29 on proprietary basis- Inviting comments thereon.

\*\*\*\*

The request has been received from Dr. Manpreet Kaur, Deptt. of Anaesthesiology,

AIIMS to purchase the subject item from M/s Vygon India Pvt. Ltd. (Mfg. by M/s Prodol

Meditech SA, Spain) on proprietary basis. The proposal submitted by M/s Vygon India Pvt.

**Ltd.** and Performa Invoice and Departmental PAC certifications are attached.

The above documents are being uploaded for open information to submit objections, comments,

if any, from any manufacturer regarding proprietary nature of the equipment/item within issue of

15 days giving reference No. 35/Prop/Anaesth/MK/2019-20/RS. The comments should be

received by office of Stores Officer (RS), Research Section at AIIMS on or before 20/01/2019

upto 12:00 p.m., failing which it will be presumed that any other vendor is having no comment

to offer and case will be decided on merits.

**STORES OFFICER (RS)** 

**Encl: Related documents enclosed.** 

1. PAC Certificate enclosed.

2. Performa Invoice



AIIMS/NR/AVP/DS/1511

Date: 15-Nov-2019

Dr. Manpreet Kaur Department of Anaesthesia AIIMS Delhi

Dear Madam,

We thank you very much for your interest in Vygon products. As desired please find our SPECIAL OFFER for the following products to meet your requirements.

Please note that the price quoted are exclusive of GST that will be charged extra at the time of delivery.

	QUOTATION		** ** ** **	GST %
S.N.	Products	HSN CODE	Unit Price (Excluding GST)	681 %
1	Airtraq Machine			
	***Technical Specification			
	A. Reusable Optics (A-501)			
	Contains the optical, anti-fog and electronic systems. The Optics element only works when it is inserted into a Blade.			
	When the battery is fully charged the Optics can be used for approximately 100 intubations.			
	Disposable Blade and Eyecup			
	· Fully protects Optics from patient contact.			
	Quick and easy assembly with Optics.			
	· Include lenses at the distal end of blade and also in eyecup.			
	Docking Station (A-590)			
	· Recharges the battery.			
	Displays the remaining number of uses of the Optics.	90189099	34132	12%
	B. Airtraq Camera (A-390K)			
	Characteristics			
	·Wireless image transmission (When used with the Airtraq Wireless			
	Display Recorder monitor)			
	Attachable to every Airtraq model			
	Connectable to any external display or monitor which has			
	an available RCA (Radio Corporation of America or Phone			
	Connector) or BNC (Bayonet Neill Concelman or Barrel Nut			
	Connector) video input jack Lightweight.	90063000	35750	18%
	C. Blades (A511/521) Total - Nos (15 M/15 F - total 50			
	blades, 25 each)	90189099	33100	12%

After the receipt for your confirmed order, the goods shall be supplied immediately or after 4 - 6 weeks in case of nonavailability The supply will be made through over authorized distributorM/S Ansh Healthcare, Delhi.

We assure you with our best services at any time.

Thanking you For Vygon India PVI. Little A PVT. LTD.

Deeapk Saluja Authorised Signatory Associate Vice President-Sales & Marketing Mob: +91 9818220820

Email Id : dsaluja@vygonindia.com

VYGON INDIA PVT. LTD.

B-17, Infocity, Sector-34, Gurgaon-122 001 (Haryana) Tel: (0124) 400 2801 / 2 Fax: (0124) 4002803 www.vygon.com GSTIN: 06AABCV0776C1ZL

Tin No.: 08391820479 IEC No.: 0500020183/07-07-2000 DL No. 12068 OW/h & 11881

Nega. Office: Vygon India Pvt. Ltd. 573, 3rd Floor, Main Road Chirag Delhi, New Delhi-110017 CIN: U74899DL1995PTC069083 PAN No.: AABCV0776C

Head Office: Head Office: Vygon. S.A. 5/11, Rue Adelfine B.P.7 - 95440 Ecouen, France Tel: (331) 3992 - 6410, Fax: (331) 3992-6444 Telex 606090





# Mowing are the specifications:

#### VIDEO LARYNGOSCOPE SPECIFICATIONS

Videolaryngoscope Machine

\*\*\*Technical Specification

#### A. VIDEOLARYNGOSCOPE OPTICS : Reusable Optics

- -Contains the optical, anti-fog and electronic systems. The Optics element only works when it is inserted into a Blade.
- ·When the battery is fully charged the Optics can be used for approximately 100 intubations.

Disposable Blade and Eyecup

- · Fully protects Optics from patient contact.
- · Quick and easy assembly with Optics.
- · Include lenses at the distal end of blade and also in eyecup.

#### **Docking Station**

- · Recharges the battery.
- ·Displays the remaining number of uses of the Optics.

#### B. VIDEOLARYNGOSCOPE WIFI CAMERA

Characteristics

- ·Wireless image transmission (When used with the Videolaryngoscope Wireless Display Recorder monitor)
- ·Attachable to every Videolaryngoscope model

Connectable to any external display or monitor which has an available RCA (Radio Corporation of America or Phone Connector) or BNC (Bayonet Neill Concelman or Barrel Nut Connector) video input jack. - Lightweight.

C. Blades

Disposables blades in two sizes 1. Videolaryngoscope Regular Blade (Size 3) ETT Size: 7.0 – 8.5 Videolaryngoscope small Blade (Size 2) ETT Size: 6.0 - 7.5

सहायक आधार्य/ASSISTANT PROFESSOR संवेदनाहरू विज्ञान संवेदनाहरू विज्ञान PARTMENT OF ANAESTHESHOLOGE जाखिल भारतीय आसुनिकान संस्थान अस्त्र अस्टारामा का महामदा अस्टारा

सहायक बाबार्व/ASSISTANT

CO-INVESTIGATORS



## PRODOL

Prodol Meditec SA Muelle T. Olavarri, 5 3rd, 48930 Las Arenas Vizcaya Spain

25th November 2017

### To Whom Soever it may Concern

Dear Sir,

This is to certify & confirm that Airtraq SP Video Laryngoscope, Airtraq Avant Video Laryngoscope, A360 WiFi camera and A-307 smart phone adapters with their unique 90 degree design and multiple optical function are being solely manufactured by Prodol Meditec SA having its manufacturing office at 1/f Block C No 18, 7th Science Ave, Hi-Tech Innovation Coast, Tangjia Bay, Zhehal, China.

This item is proprietary of Prodol Meditec SA.

For Vygon India who are our sole exclusive distributor in India

**Authorised Signatory** 

Mr. Nigel Dunlop Prodol Meditec SA

Product Manager

Nigel dunlop@airtraq.com

+44 77 987 466 96